

BIO S1

Handheld Digital Medical Scope Series



www.portable-ultrasound-scanner.com

Quality, Professional, Innovational

Handheld Digital Medical Scope Series BIO S1

Features


1. Bigger LED Screen 3.5 inch with 720*480 resolution.
2. Never hot lens with 300,000 pixels.
3. 10 pcs free disposable earmuff to prevent ear diseases infection.
4. Built-in optional wifi to ensure freely connection.
5. Examine external ear and middle ear
6. Cerumen, otitis media, ear infections, perforation of tympanic membrane checking
7. Small & Light to hold, Easy to carry with aluminum box.
8. Otoscope, Rhinoscope, Laryngoscope three lens Optional.
9. Rechargeable Li-ion battery to ensure 4 hours continuous working at least.
10. Cost-effective price












Specifications

LCD screen:	3.5 inch LCD
Light source:	Natural white LED
White balance:	Automatic color calibration
Resolution:	720 X 480
Focus range:	20~60 mm
File format:	JPG (still image) , AVI (video)
Interface:	USB,AV output
Storage:	Micro SD card
Battery:	Li-ion rechargeable battery
Power supply:	100V~240V AC50/60Hz DC IN+5V 1A

Original images took by S1



Components

 <p>S1 main unit</p>	 <p>Otoscope lens</p>
 <p>Aluminum box</p>	 <p>5V Power adaptor</p>
 <p>3.7V 18650 Li-ion battery</p>	 <p>USB charging cable</p>
 <p>AV output head</p>	 <p>TF card</p>
 <p>10 pcs Disposable ear muff</p>	

Optional:

<p>Wifi module</p>  <p>Wifi button</p>
<p>Rhinoscope</p> 
<p>Laryngoscope 5pcs disposable mouth mat</p> 
<p>Monitor</p> 

Wuxi Biomedical Technology Co., Ltd.

www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobilephone: 86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.